

LCP crystallization screen - technical sheet

Formulations: (Patent Pending)

| Well ID | <u>precipitant</u> | <u>buffer (0.1 M)</u> | <u>salt (0.2 M)</u> | Well ID |
|----------------|---------------------------|------------------------------|---------------------------------|----------------|
| A 1 | 10% (v/v) 2-propanol | Tris pH 7.0 | Zn(OAc) ₂ | A 1 |
| A 2 | 10% (v/v) 2-propanol | acetate pH 4.5 | Ca(OAc) ₂ | A 2 |
| A 3 | 10% (v/v) 2-propanol | citrate pH 5.5 | Li ₂ SO ₄ | A 3 |
| A 4 | 10% (v/v) 2-propanol | Tris pH 7.0 | Ca(OAc) ₂ | A 4 |
| A 5 | 10% (v/v) 2-propanol | acetate pH 4.5 | NaCl | A 5 |
| A 6 | 10% (v/v) 2-propanol | citrate pH 5.5 | MgCl ₂ | A 6 |
| A 7 | 10% (w/v) PEG-8000 | HEPES pH 7.5 | NaCl | A 7 |
| A 8 | 10% (w/v) PEG-8000 | MES pH 6.0 | Zn(OAc) ₂ | A 8 |
| A 9 | 10% (w/v) PEG-8000 | Na/K phosphate pH 6.2 | Li ₂ SO ₄ | A 9 |
| A 10 | 10% (w/v) PEG-8000 | acetate pH 4.5 | Ca(OAc) ₂ | A 10 |
| A 11 | 10% (w/v) PEG-8000 | Tris pH 8.5 | MgCl ₂ | A 11 |
| A 12 | 10% (w/v) PEG-8000 | cacodylate pH 6.5 | NaCl | A 12 |
| B 1 | 10% (w/v) PEG-8000 | citrate pH 5.5 | Li ₂ SO ₄ | B 1 |
| B 2 | 10% (w/v) PEG-8000 | HEPES pH 7.5 | MgCl ₂ | B 2 |
| B 3 | 20% (w/v) PEG-2000 MME | HEPES pH 7.5 | NaCl | B 3 |
| B 4 | 20% (w/v) PEG-2000 MME | MES pH 6.0 | Zn(OAc) ₂ | B 4 |
| B 5 | 20% (w/v) PEG-2000 MME | Na/K phosphate pH 6.2 | Li ₂ SO ₄ | B 5 |
| B 6 | 20% (w/v) PEG-2000 MME | acetate pH 4.5 | Ca(OAc) ₂ | B 6 |
| B 7 | 20% (w/v) PEG-2000 MME | Tris pH 8.5 | MgCl ₂ | B 7 |
| B 8 | 20% (w/v) PEG-2000 MME | cacodylate pH 6.5 | Li ₂ SO ₄ | B 8 |
| B 9 | 20% (w/v) PEG-2000 MME | HEPES pH 7.5 | Li ₂ SO ₄ | B 9 |
| B 10 | 20% 1,4-butanediol | Tris pH 8.5 | Li ₂ SO ₄ | B 10 |
| B 11 | 20% 1,4-butanediol | cacodylate pH 6.5 | Ca(OAc) ₂ | B 11 |
| B 12 | 20% 1,4-butanediol | citrate pH 5.5 | MgCl ₂ | B 12 |
| C 1 | 20% 1,4-butanediol | Tris pH 7.0 | Zn(OAc) ₂ | C 1 |
| C 2 | 20% 1,4-butanediol | Tris pH 8.5 | MgCl ₂ | C 2 |
| C 3 | 20% 1,4-butanediol | citrate pH 5.5 | NaCl | C 3 |
| C 4 | 20% (w/v) PEG-1000 | HEPES pH 7.5 | NaCl | C 4 |
| C 5 | 20% (w/v) PEG-1000 | MES pH 6.0 | Zn(OAc) ₂ | C 5 |
| C 6 | 20% (w/v) PEG-1000 | Tris pH 7.0 | Li ₂ SO ₄ | C 6 |
| C 7 | 20% (w/v) PEG-1000 | HEPES pH 7.5 | MgCl ₂ | C 7 |
| C 8 | 20% (w/v) PEG-1000 | MES pH 6.0 | NaCl | C 8 |
| C 9 | 2.5 M NaCl | HEPES pH 7.5 | | C 9 |
| C 10 | 2.5 M NaCl | MES pH 6.0 | Zn(OAc) ₂ | C 10 |
| C 11 | 2.5 M NaCl | Tris pH 8.5 | Li ₂ SO ₄ | C 11 |
| C 12 | 2.5 M NaCl | cacodylate pH 6.5 | Ca(OAc) ₂ | C 12 |
| D 1 | 2.5 M NaCl | citrate pH 5.5 | MgCl ₂ | D 1 |
| D 2 | 2.5 M NaCl | cacodylate pH 6.5 | Li ₂ SO ₄ | D 2 |

| Well ID | <u>precipitant</u> | <u>buffer (0.1 M)</u> | <u>salt (0.2 M)</u> | Well ID |
|----------------|--|------------------------------|---------------------------------|----------------|
| D 3 | 30% (w/v) PEG-400 | Tris pH 7.0 | Zn(OAc) ₂ | D 3 |
| D 4 | 30% (w/v) PEG-400 | Imidazole pH 8.0 | Li ₂ SO ₄ | D 4 |
| D 5 | 30% (w/v) PEG-400 | citrate pH 5.5 | MgCl ₂ | D 5 |
| D 6 | 30% (w/v) PEG-400 | citrate pH 5.5 | Li ₂ SO ₄ | D 6 |
| D 7 | 30% (w/v) PEG-400 | Imidazole pH 8.0 | Zn(OAc) ₂ | D 7 |
| D 8 | 30% (w/v) PEG-400 | Tris pH 7.0 | NaCl | D 8 |
| D 9 | 15% (v/v) Ethanol | Na/K phosphate pH 6.2 | NaCl | D 9 |
| D 10 | 15% (v/v) Ethanol | acetate pH 4.5 | Zn(OAc) ₂ | D 10 |
| D 11 | 15% (v/v) Ethanol | cacodylate pH 6.5 | Ca(OAc) ₂ | D 11 |
| D 12 | 15% (v/v) Ethanol | acetate pH 4.5 | MgCl ₂ | D 12 |
| E 1 | 15% (v/v) Ethanol | cacodylate pH 6.5 | NaCl | E 1 |
| E 2 | 10% (w/v) PEG-3000 | HEPES pH 7.5 | NaCl | E 2 |
| E 3 | 10% (w/v) PEG-3000 | MES pH 6.0 | Zn(OAc) ₂ | E 3 |
| E 4 | 10% (w/v) PEG-3000 | Tris pH 8.5 | Ca(OAc) ₂ | E 4 |
| E 5 | 10% (w/v) PEG-3000 | citrate pH 5.5 | MgCl ₂ | E 5 |
| E 6 | 10% (w/v) PEG-3000 | Tris pH 8.5 | Li ₂ SO ₄ | E 6 |
| E 7 | 10% (w/v) PEG-3000 | MES pH 6.0 | Li ₂ SO ₄ | E 7 |
| E 8 | 1.0 M (NH ₄) ₂ HPO ₄ | HEPES pH 7.5 | NaCl | E 8 |
| E 9 | 1.0 M (NH ₄) ₂ HPO ₄ | cacodylate pH 6.5 | NaCl | E 9 |
| E 10 | 1.26 M (NH ₄) ₂ SO ₄ | Na/K phosphate pH 6.2 | Li ₂ SO ₄ | E 10 |
| E 11 | 1.26 M (NH ₄) ₂ SO ₄ | Tris pH 7.0 | MgCl ₂ | E 11 |
| E 12 | 1.26 M (NH ₄) ₂ SO ₄ | Imidazole pH 8.0 | Li ₂ SO ₄ | E 12 |
| F 1 | 1.26 M (NH ₄) ₂ SO ₄ | citrate pH 5.5 | NaCl | F 1 |
| F 2 | 20% (w/v) PEG-8000 | HEPES pH 7.5 | NaCl | F 2 |
| F 3 | 20% (w/v) PEG-8000 | Tris pH 7.0 | Zn(OAc) ₂ | F 3 |
| F 4 | 20% (w/v) PEG-8000 | Imidazole pH 8.0 | Li ₂ SO ₄ | F 4 |
| F 5 | 20% (w/v) PEG-8000 | acetate pH 4.5 | MgCl ₂ | F 5 |
| F 6 | 20% (w/v) PEG-8000 | cacodylate pH 6.5 | NaCl | F 6 |
| F 7 | 20% (w/v) PEG-8000 | Tris pH 7.0 | MgCl ₂ | F 7 |
| F 8 | 1.0 M Sodium citrate | HEPES pH 7.5 | NaCl | F 8 |
| F 9 | 1.0 M Sodium citrate | MES pH 6.0 | Zn(OAc) ₂ | F 9 |
| F 10 | 1.0 M Sodium citrate | Imidazole pH 8.0 | Li ₂ SO ₄ | F 10 |
| F 11 | 1.0 M Sodium citrate | acetate pH 4.5 | MgCl ₂ | F 11 |
| F 12 | 1.0 M Sodium citrate | Tris pH 8.5 | Li ₂ SO ₄ | F 12 |
| G 1 | 1.0 M Sodium citrate | MES pH 6.0 | MgCl ₂ | G 1 |
| G 2 | 10% (v/v) 2-propanol | acetate pH 4.5 | Li ₂ SO ₄ | G 2 |
| G 3 | 10% (v/v) 2-propanol | citrate pH 5.5 | NaCl | G 3 |
| G 4 | 10% (w/v) PEG-8000 | acetate pH 4.5 | MgCl ₂ | G 4 |
| G 5 | 10% (w/v) PEG-8000 | Tris pH 8.5 | NaCl | G 5 |
| G 6 | 20% (w/v) PEG-2000 MME | acetate pH 4.5 | Zn(OAc) ₂ | G 6 |
| G 7 | 20% (w/v) PEG-2000 MME | Tris pH 8.5 | Li ₂ SO ₄ | G 7 |
| G 8 | 20% 1,4-butanediol | Tris pH 8.5 | NaCl | G 8 |
| G 9 | 20% 1,4-butanediol | cacodylate pH 6.5 | Li ₂ SO ₄ | G 9 |
| G 10 | 20% (w/v) PEG-1000 | HEPES pH 7.5 | Li ₂ SO ₄ | G 10 |

| Well ID | <u>precipitant</u> | <u>buffer (0.1 M)</u> | <u>salt (0.2 M)</u> | Well ID |
|----------------|---------------------------|------------------------------|---------------------------------|----------------|
| G 11 | 20% (w/v) PEG-1000 | MES pH 6.0 | MgCl ₂ | G 11 |
| G 12 | 2.5 M NaCl | MES pH 6.0 | Li ₂ SO ₄ | G 12 |
| H 1 | 30% (w/v) PEG-8000 | Tris pH 8.5 | Li ₂ SO ₄ | H 1 |
| H 2 | 30% (w/v) PEG-8000 | cacodylate pH 6.5 | Ca(Oac) ₂ | H 2 |
| H 3 | 30% (w/v) PEG-400 | MES pH 6.0 | MgCl ₂ | H 3 |
| H 4 | 30% (w/v) PEG-400 | Tris pH 7.0 | Li ₂ SO ₄ | H 4 |
| H 5 | 10% (w/v) PEG-3000 | HEPES pH 7.5 | Li ₂ SO ₄ | H 5 |
| H 6 | 10% (w/v) PEG-3000 | MES pH 6.0 | MgCl ₂ | H 6 |
| H 7 | 1.0 M Sodium citrate | Imidazole pH 8.0 | MgCl ₂ | H 7 |
| H 8 | 1.0 M Sodium citrate | Tris pH 8.5 | NaCl | H 8 |
| H 9 | 2.5M Sodium Malonate | Tris pH 7.0 | | H 9 |
| H 10 | 2.5M Sodium Malonate | acetate pH 4.5 | | H 10 |
| H 11 | 2.5M Sodium Malonate | Tris pH 8.5 | | H 11 |
| H 12 | 2.5M Sodium Malonate | Imidazole pH 8.0 | | H 12 |