



Emerald BioSystems

Cubic 96 WELL BLOCK

The format for the 96 block well plates is as follows:

A1 through A12 contains Cubic conditions **1-12**

B1 through B12 contains Cubic conditions **13-24**

C1 through C12 contains Cubic conditions **25-36**

D1 through D12 contains Cubic conditions **37-48**

E1 through E12 contains Cubic conditions **49-60**

F1 through F12 contains Cubic conditions **61-72**

G1 through G12 contains Cubic conditions **73-84**

H1 through H12 contains Cubic conditions **85-96**