


Full Ozma PEG-Series in Block Plates



Catalogue No.: EBS-PEG-48BLK Full Ozma™ PEG-Series 48-Salt in 2 x 96 Well Matrix Block Plates: Contains 192 unique formulations comprising all 4 sets of our Ozma PEG 1K, 4K, 8K, and 10K screens, each with 48 different salts.

Used for the crystallization of biological macromolecules in PEG-1000, PEG-4000, PEG-8000 and PEG-10,000 by varying the salt type. 1.0 ml each prepared with ultra-pure chemicals and water (18.2 Megohm-cm) followed by sterile 0.22 micron filtration into a sterile 96 deep well matrix block plate that is sealed with a sterile reusable Matte cover.

Use in combination buffers covering pH ranges from 2.4 to 11.5 that comprise our pH at Buffer Screen to further expand the crystallization possibilities.